



FEATURES

1. Narrow & Speech band masking also available with wide masking.
2. Supplied with HD-01 as standard, can accept both 8 and 10 ohms headphones. Optional TDH 39/49 and DD 45 available at extra cost
3. TONE DECAY test available along with pulse tone & +10dB facility.
4. Digital Calibration.
(Can be calibrated from front panel without moving any part inside the audiometer.)
5. With USB port for computer adaptability to operate through Computer Software provided Free of charge.
6. Warranty of 24 months. (Accessories not included)
7. Free Field option included.
8. AC Left & Right, BC Left & Right, Free field, Masking, Calibration done separately.

SPECIFICATIONS

Input Voltage: AC 230 Volts - 50 Hz/DC 12 Volts (optional)
 Pure Tone Frequencies:
 Air Conduction : 125 Hz to 8 KHz (Maximum 120 dB)
 Speech : 100dB
 Bone Conduction: 250Hz to 8KHz (Maximum 70 dB)
 Speech: 70dB
 Noise: Wide, Narrow & Speech Band
 Attenuator (Masking) In steps of 5dB
 Masking Range 0 dB to 100 dB.
 Weight: 1.5 kg Nett.. With bag 4 kg (Approx)
 Dimension in cm : 20 cm x 27 cm x 8.5 cm
 (Specs- Subject to change due to continuous development)
 Complaint with ANSI S3.6 1989 standard
 (if used with Standard Accessories)



GEM SURGITECH
 MANUFACTURER OF MEDICAL EQUIPMENTS & INSTRUMENTS

Factory: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India
 Contact Info: Ph No. +91-9996256235, +91-9416025292
 Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com
 Website: www.gemsurgitech.com
 (An ISO 9001:2015 Certified Company)