



- ☐ 4.5 mm Ø
- ☐ 5.5 mm Ø
- ☐ 6.5 mm Ø
- ☐ 7.5 mm Ø

O.001
GRUBER
Aural Specula
Set of 4



- ☐ 4.5 mm Ø
- ☐ 5.5 mm Ø
- ☐ 6.5 mm Ø
- ☐ 7.5 mm Ø

O.002
HARTMANN
Aural Specula
Set of 4



O.003
HOLMGREN
Aural Specula



O.004
ROSEN
Aural Specula
Set of 4



O.005
SHEA
Aural Specula
Set of 4



O.006
SIEGLE
Pneumatic Speculum
Set of 4



O.007ROSEN
Aural Slotted Specula
Set of 4



O.008
TUMARKIN
Tapered Slotted Aural
Specula, Set of 4



O.009
Aural Syringe All Metal
with 2 Pipes Conical and
Bulbous Pipes



O.010
HARTMANN
Tuning Forks
128-256-512-1024-
2048-4096



O.011
GARDINER
Tuning Forks
128-256-512-1024-
2048-4096



O.012
EUSTACHIAN
Catheters




O.013
FORMBY
Cerumen Scoop
and Hook



O.014
JOBSON
Horn Probes



O.015
HOVELL
Loop



O.016
Ear Victus
with Scoop



O.016
Foreign Body
Removing



O.018
Foreign Body
And Hard Was
Removing



O.019
WILDE
Aural Forceps



O.020
FAGGE
Aural Forceps



O.021
ADSON
Plain Forceps




O.022
ADSON
Forceps with
1 x 2 teeth



O.023
OSSICLE
Holding
Forceps



O.024
HARTMANN
Aural Forceps



O.025
TILLEY
Aural
Forceps



O.026
LEMPERT
Endaural
Speculum



O.027
LEMPERT Rongeur,
Straight



O.028
LEMPERT Rongeur,
Curved



O.029
LEMPERT Rongeur,
Double Action



O.030
PLASTER Self
Retaining Retractor
2 x 2 Prongs



O.031
WULLSTEIN
Self Retaining Retractor
3 x 3 Prongs



O.032
WULLSTEIN
Self Retaining Retractor
3 x 4 Prongs



O.033
Self Retaining
Adjustable Retractor



O.034
Jansen
Retractor



O.035
Plaster Jansen Retractor
Left Solid Blade



O.035B
Plaster Jansen Retractor
Right Solid Blade



O.036
Lempert Mastoid
Curette



O.037
GLEGG Aural Snare



O.038
Jansen Chisel



O.039
Jansen Osteotome
2,4,6,8,10mm



O.040
Jenkins GOUGE,
2,4,6,8,10mm



O.041
FARABEUF
Elevator Straight



O.041B
FARABEUF
Elevator Curved



O.042
Kerrison Punch
Upward Cutting



O.043
Lake Seeker



O.044
Houses Graft
Press Forceps



O.045
Lempert Mastoid Suction Tube
Size 1-2-3-4



O.046
Ferguson Suction
Tube Size 1-2-3-4



O.047
Verhoeven Suction Tube SWG 14,
15, 16 17, 18, 20, 21, 22, 23, 24



O.048
VERHOEVEN
Suction Tube
Holder



O.049
Tube Connector




O.050
Temporal Bone Holder



Micro Aural Crocodile
Forceps Serrated

- O.051 Fine Straight
- O.052 Fine 1 & 2 Teeth
- O.053 Right Curved
- O.054 Left Curved



Micro Aural Cup Forceps

- O.055 Straight Cup
- O.056 Right Curved
- O.057 Left Curved



Micro Aural Scissors

- O.058 Right Curved
- O.059 Left Curved
- O.060 Straight Fine



O.061
Mc-Gee Wire Closing
Forceps



O.062
Teflon Piston Holding
Forceps



O.063
Grommet Inserter



O.064
Piston Introducing Forceps



O.065 Up
O.066 Down



O.067

O.067 WULLSTEIN Needle Straight, Sharp

O.068 WULLSTEIN Needle Slight Curved, Sharp

O.069 WULLSTEIN Needle Strong Curved, Sharp

O.070 WULLSTEIN Needle Very Slight Curved, Sharp

O.071 MICRO PICK 90° Right Angled

O.072 MICRO PICK 45° Angled

O.073 PLESTER Foot Plate Hook

O.068

O.069

O.070

O.071

O.072

O.073

O.074

O.074 SCHUHKNECHT Micro Needle Straight, Sharp

O.075 SCHUHKNECHT Micro Needle 90° Right Angled, Right

O.076 SCHUHKNECHT Micro Needle 90° Right Angled, Left

O.077 SCHUHKNECHT Micro Needle 45° Right Angled, Left

O.078 SCHUHKNECHT Micro Needle 45° Right Angled, Right

O.075

O.076

O.077

O.078

O.079
ROSEN Elevator
Wide LengthO.080
BEALES
RaspatoryO.081
Sickle KnifeO.082
PLESTER
First Incision KnifeO.083
MYRINGOTOMY
Knife



O.084
Rosen Curette



O.085
Rosen Curette
Double Curved



O.086
Rosen Curette
Curved



O.087
Circular Cutting
Knife 60°



O.088
Circular Cutting
Knife Curved



O.089
House Curette
Double Ended



O.090
LARKIN Hand
Trephine



O.091
LARKIN
Penetrator (Burr)



O.092
LARKIN
Perforator



O.093
House Measuring
Roads 3.5, 4.0,
4.5mm




O.094
Hook Ball
End Bent 45°



O.095
SHEA Piston
Depth Gauge
3.5, 4.0, 4.5mm




O.096
Fine Hook




O.097
Fine Knife
Downward
Cutting




O.098
Fine Knife
Upward
Cutting



O.099
Spear
Pointed
Fine Knife
Cutting Left




O.100
Spear
Pointed
Fine Knife
Cutting Right




O.101
Raspatory
Fine
Curved Left




O.102
Raspatory
Fine
Curved Right




O.103
Cutting Burr
SS



7.0



6.0




5.0



4.0




3.1



2.7



1.8



O.104
Tungstone
Carbide
Cutting Burr



7.0



6.0



5.0



4.0



3.1



2.7



1.8




O.105
Diamond
Burr



7.0



6.0



5.0



4.0



3.1



2.7



1.8



O.106
Oval Shape
Burr SS



7.0



6.0



5.0



BURR STAND BLOCK



O.108
Hand Piece
Straight



O.109
Hand Piece
Angled



O.110
Hand Piece
Spray



O.111
Hand Piece for
Wrench

**ENT SURGICAL MICROSCOPE (Make:GEM)**

Serial No. O.113

Technical Specifications:-

- Head: Straight Binocular Tube
- Magnification: 0.4x, 0.6x, 1.0x, 1.6x, 2.5x
- Eye Pieces: SWF 12.5x paired super wide field
- Working Dist: F=200mm Objective Lens
- Fine Focusing: Manual Fine Focusing
- Light Source: 100W LED Light Source with bright white LED Light
- Arm: Counter balanced pantographic arm with 320° rotation
- Floor Stand: Mobile floor stand on four caster wheels for easy handling and absolute stability
- Power Supply: AC 220V – 240V

Optional Accessories:-

- Tilttable/Inclined Head
- Beam Splitter with 'C' Mount
- CCD Digital Camera
- LED Screen for Live Surgery
- Motorised with Foot Control
- Objective lens: 250mm, 300mm & 400mm



O.117
Sterilization and Storage Micro
Instrument Case



O.118 Binocular Loupe



O.120 Micro Instrument Case
with Silicone Mat



O.121 Rack for holding Aural
Micro Instruments



O.122 Micro Instruments Box



R.001
HARTMANN
Nasal Speculum



R.002
KILLIAN Nasal
Speculum 5 cm,
6.5cm, 7.5 mm, 9cm



R.003
PILCHER Nasal
Speculum



R.004
KILLIAN Nasal Speculum with
Fibre Optic Light carrier



R.005
THUDICHUM
Nasal Speculum, Sec of 3



R.006
COTTLE
Retractor Blunt



R.007
COTTLE
Retractor
Sharp



R.008
GILLIES
Skin Hook



R-009
COTTLE Uper
Lateral Exposing
Retractor



R-010
ALAR
Retractor



R-011
KILNER
Retractor



R-013
COTTLE Elevator
Slightly Curved



R-014
JOSEPH Nasal
Elevator



R-015
FREER Elevator
Double Ended



R-016
Freer
Elevator



R-017
MASING Elevator
With Graduated
Freer Elevator



R-018
SMR Elevator
One Side Straight
One Side Curved



R-019
SMR Fine
Elevator



R-020
HOWARTH
Elevator



R-021
PARKER
Elevator



R-022
KILLIAN
Nasal
Gauge



R-023
TILLEY Septum
Gauge V Shape
Cutting Edge



R-024
BALLENGER Swivel
Knife Bayonet



R-025
KRAUSE Nasal
Snare



R-026
HEYMANN
Turbinectomy Scissors



R-027
LUC Nasal Turbinate
Forceps Oval



R-028
ZOELLNER
Suction Tube



R-029
LUC Nasal Turbinate
Forceps, Small Size



R-030
COTTLE Septal
Scissors Spring
Angular Bows



R-031
FOMAN Angular
Scissors



R-032
LEMPERT
Suction and
Irrigation Tubes
Size 1-2-3-4



R-033
FERGUSON
Suction Tube,
Size 1-2-3-4-5



R-034
KNAPP Scissor
Straight



R-035
KNAPP Scissor
Curved



R-036
Fine Dressing
Scissor Straight



R-037
Fine Dressing
Scissors Curved



R-038
FOMAN Upper
Lateral Scissor



R-039
TILLEY LICHTWITZ
Antral Trocar



R-040
Skin Biopsy Punch
Forceps



R-041
ST. BARTHOLOMEWS
Ethmoid Punches
Straight



R-042
ST. BARTHOLOMEWS
Ethmoid Punches
Upward



R-043
Antral Wash Rubber
Syringe



ABBEY Sub Mucosal Resection Needle,
Knife, Suction Tube

R-044 Knife

R-045 Fine Needle

R-046 Fine Suction Tube with Handle





R-047
WATSON WILLIAMS
Nasal Polypus Forceps



R-048
WATSON WILLIAMS
Forceps Fine Serrated
15 cm



R-049
WAUGH Dissecting
Forceps Fine 1 x2
Teeth, 15 cm



R-050
Septum Cutting
Forceps



R-051
WALSHAM Septum
Straightening Forceps



R-052
WALSHAM Septum
Straightening Forceps
Left **R-053** Right



R-054
ASH Septum
Straightening Forceps



R-055
Septum Morcelizer
Double Action




R-056
TILLEY Harpoon
Trocar



R-057
MYLES Neso
Antral Perforator



R-058
TILLEY Antrum
Burs



R-059
WIENER Antrum
Rasp



R-060
KERRISON Punch
Upward Cutting



R-061
AUFRIGHT
Nasal Retractor



R-062
Silver Chisel
Right



R-063
Silver Chisel
Left



R-064
Silver Chisel
Straight



R-065
BOZNER Nasal
Osteotome 6mm,
8mm, 10mm, 12mm



R-066
Mallet/
Hammer



R-067
GILLIES
Chisel 2mm,
4mm, 6mm,
8mm, 10mm,



R-067
GILLIES
Osteotome
2mm, 4mm, 6mm,
8mm, 10mm,



R-067
GILLIES Gouge
2mm, 4mm,
6mm, 8mm,
10mm,



R-070
MALTZ
Rasp



R-071
AUFRICHT
Nasal Rasp



R-072
JOSEPH
Nasal
Rasp



R-073
JOSEPH
Nasal
Rasp Fine



R-074
JOSEPH
Nasal
Saw Left




R-075
JOSEPH
Nasal Saw
Right



R-076
Caliper



R-077
KNAPP Columnella
Forceps




R-078
Wire Cutting
Scissors




R-079
HUMBY Skin
Graft Knife




R-080
Silver Skin Grafting
Knife For Saving
Blades Use

**S-001**


Sinsucope 4mm 0°, wide angle,
Working length 175mm

**S.002**


Sinsucope 4mm 30°, wide angle,
Working length 175mm

**S.003**

Sinsucope 4mm 45°, wide angle,
Working length 175mm

**S.003A**

Sinsucope 4mm 70°, wide angle,
Working length 175mm

**S.004**

Sinuscopy Handle
Round/Flat 11cm
length



Trocar and Cannula
Sinuscopy,

S.005 Billshaped
beak

S.006 Oblique beak



VELCKEN
Antrum
Cannula, 15
cm,
S.007



Bipolar
Suction
Tube 15 cm,
S.008



FRAZIER
Suction Tube,
ball ended
S.009



Bipolar
Suction Tube,
ball ended
S.010



Elevator with
Knife
S.011



S.012
Freer
Elevator



S.013
Masing
Elevator



S.014
Antrum
Curette



S.015
Antrum
Curette



S.016
Antrum
Ball probe



S.017
Antrum Ball
probe Double
Ended



S.018
Antrum
Curette



S.019
Elevator/Curette
Double Ended



S.020
Sickle Knife



S.021
Straight



S.022
45° UP



S.023
90° UP

Blakesley Nasal Forceps




S.024
UP



S.025
Down

Antrum Punch Backward Cutting Forceps



S.026
Straight



S.027
Right



S.028
Left

Nasal Scissors



S.029
Straight



S.030
45° UP

Through Cutting Punch Forceps





S.031



S.032




Frontal Sinus Punch Forceps (Giraff) 55° Double Action



S.033



Through Cutting Punch Forceps Upward (Rotatable)

S.034
1mmS.035
2mmS.036
3mm

Antrum Punch Through Cutting 13cm

S.037
1mmS.038
2mmS.039
3mm

D.C.R. Punch



T.001
Tongue
Depressor

T.002

T.003



T.004
Birkett Tonsil Artery
Forceps Straight



T.005
Birkett Tonsil Artery
Forceps Curved



T.006
Negus Tonsil Artery
Forceps



T.007
Negus Tonsil Artery
Forceps



T.008
Woods Tonsil
Scissors Straight



T.009
Woods Tonsil
Scissors Curved



T.010
Denis Browne Tonsil
Holding Forceps



T.011
Colver Tonsil Holding
Forceps



T.012
Tonsil Needle Holder



T.013
Peritonsillar
Forceps



T.014
EVE Tonsil
Snare



T.015
Irwin Moore
Tonsil Needle
Right & Left



T.016
Tonsillar Knire



T.017
Hanke Tonsil
Dissector -
One Side
Serreted



T.018
Hurd Tonsil
Dissector & Pillar
Retractor



T.019
Tonsil Pillar
Retractor
& Dissector
Serrated Fine



T.020
Negus
Pusher
Knot Tier



T.021
Tongue
Tie
Probe



T.022
Yankauer
Tonsil
Scution
Tube



T.023
Heister Mouth
Gag



T.024
Jennings
Mouth Gag



T.025
Waugh Tonsil
Dissection
Forceps



T.026
Waugh Tonsil
Dissection
Forceps 1 x 2
Teeth



T.027
Backmann
Adenoid
Curette



T.028
S.T Clair Thompson
Adenoid Curette
with cage



T.029
Davis Boyle Mouth
Gag with Tonsil
Plates



T.030
Doughty Slotted Mouth
Gag with Tongue Plates



T.031
Dingman Mouth Gag
with Tongue Plates



T.032
Draffin Bipod Stand
with Plate