

ENT LED HEADLIGHT

MODELS: GMH3/5/7



GEM SURGITECH
MANUFACTURER OF MEDICAL EQUIPMENTS & INSTRUMENTS



Main Features:

- Lightweight
- Soft Foam Inserts for Comfortable Long Duration Use
- Bright and Energy Efficient LED
- Adjustable Brightness
- 2-way Adjustable Headband Perfectly Fits on the Head

Specifications:

Models:	GMH3	GMH5	GMH7
Power:	5W	5W	5W
Battery capacity:	2600mah	2600mah	2600mah
Working duration:	4hrs.	4hrs.	4hrs.
Color temperature:	6000k	6000k	6000k
Adjustable Spot Light:	N/A	N/A	Yes
Interchangeable Battery Box:	N/A	Yes	Yes
Lux:	>100000lux	>100000lux	>100000lux
Weight:	200g.	200g.	210g.

Factory: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India
Contact Info: Ph No. +91-9996256235, +91-9416025292
Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com
Website: www.gemsurgitech.com
(An ISO 9001:2015 Certified Company)