




Standard Features

- Built-in Twin Halogen Light Source 150w with Brightness Control.
- Built-in Suction Machine with Pressure Gauge (Twin Jars)
- Slim LED X-Ray View Box.
- Thermal Cautery.
- Anti Fog Device/Mirror Warmer.
- LED Pen Light for Endoscope Illumination
- X2 Stainless Steel Instruments Trays with Cover.
- X2 Medicine Bottles.
- X2 Sliding Drawers for Storage
- Head Light Band Hanger.
- X2 Endoscope Holders.
- Otoscope Holder.
- Waste Container/Cotton Gauze Box.



Colors Available:

Blue, Black, Gray, Orange

(any other color available on order subject to availability)

Optional Features

Compressed Air System (Spray System)

- X2 Spray Guns – One Straight and One Curved.
- X2 Pressure Gauges for Both Sprays.
- Convenient Single Handed Operation.

Advanced Ear Irrigation System

- Electronic Temperature Controller for Adjusting Water Temperature.
- Adjustable Pressure of Syringe Water.
- Built-in Water Tank.
- Foot Controlled for Easy Operation.

LED Light Source – 120w (inbuilt)

- Extremely Powerful 120w LED Light Source
- With Brightness Control.
- Produces upto 5,000 lumens.
- Color Temperature of 5700k-6500k
- Reliable LED Lifespan of more than 60,000 Hours.
It will replace one halogen light source outlet*

LCD Monitor with Stand

- Full HD Resolution Monitor 18"
- Adjustable Stand for Different Viewing Angles