



Technical Data:

Head Tube:	Straight Binocular Tube Optional: Ergonomic Tilttable Binocular Tube 45deg. Inclined Binocular Tube
Magnification:	5-Step Magnification: 0.4x, 0.6x, 1.0x, 1.6x, 2.5x Optional: 3-Step Magnification Step-Less Zoom Magnification
Eye Pieces:	SWF 12.5x, Super Wide Field
Working Distance:	F=200mm Objective Lens
Fine Focusing:	Manual Fine Focusing
Illumination:	100W LED Light Source with Super Bright White LED as Standard
Arm:	Counter Balance Pantographic Arm with 320deg. Rotation
Floor Stand:	Mobile Floor Stand on Four Caster Wheels for Easy Handling and Absolute Stability
Power Supply:	AC 220V-240V

Optional Accessories:

- Beam Splitter with 'C' Mount
- CCD Digital Camera
- LED Screen for Live Surgery
- Motorised with Foot Control
- Objective lens: 250mm, 300mm & 400mm