

ENT TREATMENT UNIT

MODEL: ETU-7000



GEM SURGITECH
MANUFACTURER OF MEDICAL EQUIPMENTS & INSTRUMENTS



Main Features:

- High Quality 304 Stainless Steel Body, Non Rusting, Non Magnetic Mounted on Caster Wheels
- Auto start/stop Suction with Gauge and Twin Jars
- Auto start/stop Compressed Air System - 2 Sprays for liquid Medicine, 1 Straight and 1 Curved Tip
- Auto start/stop Ear Irrigation System with Temperature Control Range Ambient to 100deg Celsius
- Built-in 100w LED Light Source with Brightness Control
- Mirror Preheat system
- LED X-Ray view Box
- Instrument Tray with Dust Cover
- Stainless Steel Gauze Dispenser
- Writing Desk
- Otoscope/Endoscopes Holder
- Headlight Hanger
- Vibration Free Adjustable Stand for Monitor
- Cotton/Gauze Container

Optional:

- ENT Microscope
- Medical Grade Monitor
- Patient Chair with Doctor's Stool
- ENT Endoscopes
- Endoscopy Camera