

LED LIGHT SOURCE

MODEL: GM-120



GEM SURGITECH

MANUFACTURER OF MEDICAL EQUIPMENTS & INSTRUMENTS



Main Features:

- Extremely Bright Medical Speciality LED (USA Imported LED)
- Highly Reliable, It has the Longest Life Span Compared to Xenon or Halogen Lamps.
- Completely Replaces 350W Xenon Light Source
- Lightweight

Specifications:

LED Power:	120W
Life Span:	Upto 60,000 hrs.
Connection Type:	Storz Fitting
Color Temperature:	6500k (White)
Weight:	2.5kg.

Factory: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India
Contact Info: Ph No. +91-9996256235, +91-9416025292
Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com
Website: www.gemsurgitech.com
(An ISO 9001:2015 Certified Company)