



Main Features:

- Powerful and Excellent Performance with Low Noise & Vibration.
- 4 Stage Speed Control
- Forward & Reverse Direction Switch
- Foot Control for Start/Stop

Specifications:

Electric Power:	40W
Speed:	Max 35,000 RPM
Dimensions:	W112xD147xH77 (mm)
Weight:	1.16kg

Contents: Control Unit, Hand Piece, Foot Control, Hand Piece Stand

Hand Piece Options:

- Straight Hand Piece
- Angled Hand Piece
- Reciprocating Hand Piece
- Sagittal Hand Piece