



Application:

EndoView *Lite* is a High Resolution and Economical Endoscopy Camera used in Basic ENT OPD Endoscopy, Student Training Programs & Other Imaging Applications.

Features:

Sony CCD Image Sensor
Lightweight and Compact
Anti-Glare Coated Sensor

Specifications:

Image Sensor:	1/3" Sony CCD Image Sensor
Resolution (H):	800TVL
Signal:	PAL (NTSC Optional)
White Balance:	Preset
Backlight Compensation:	Auto
Power Supply:	DC12V
Weight:	Approx. 85g