



- 4.5 mm Ø
- 5.5 mm Ø
- 6.5 mm Ø
- 7.5 mm Ø

O.001
GRUBER
Aural Specula
Set of 4



- 4.5 mm Ø
- 5.5 mm Ø
- 6.5 mm Ø
- 7.5 mm Ø

O.002
HARTMANN
Aural Specula
Set of 4



O.003
HOLMGREN
Aural Specula



O.004
ROSEN
Aural Specula
Set of 4



O.005
SHEA
Aural Specula
Set of 4



O.006
SIEGLE
Pneumatic Speculum
Set of 4



O.007ROSEN
Aural Slotted Specula
Set of 4



O.008
TUMARKIN
Tapered Slotted Aural
Specula, Set of 4



O.009
Aural Syringe All Metal
with 2 Pipes Conical and
Bulbous Pipes



O.010
HARTMANN
Tuning Forks
128-256-512-1024-
2048-4096



O.011
GARDINER
Tuning Forks
128-256-512-1024-
2048-4096



O.012
EUSTACHIAN
Catheters




O.013
FORMBY
Cerumen Scoop
and Hook



O.014
JOBSON
Horn Probes



O.015
HOVELL
Loop



O.016
Ear Victor
with Scoop



O.016
Foreign Body
Removing



O.018
Foreign Body
And Hard Was
Removing



O.019
WILDE
Aural Forceps



O.020
FAGGE
Aural Forceps



O.021
ADSON
Plain Forceps



O.022
ADSON
Forceps with
1 x 2 teeth



O.023
OSSICLE
Holding
Forceps



O.024
HARTMANN
Aural Forceps



O.025
TILLEY
Aural
Forceps



O.026
LEMPERT
Endaural
Speculum



O.027
LEMPERT Rongeur,
Straight



O.028
LEMPERT Rongeur,
Curved



O.029
LEMPERT Rongeur,
Double Action



O.030
PLASTER Self
Retaining Retractor
2 x 2 Prongs



O.031
WULLSTEIN
Self Retaining Retractor
3 x 3 Prongs



O.032
WULLSTEIN
Self Retaining Retractor
3 x 4 Prongs



O.033
Self Retaining
Adjustable Retractor



O.034
Jansen
Retractor



O.035
Plaster Jansen Retractor
Left Solid Blade



O.035B
Plaster Jansen Retractor
Right Solid Blade



O.036
Lempert Mastoid
Curette



O.037
GLEGG Aural Snare



O.038
Jansen Chisel



O.039
Jansen Osteotome
2,4,6,8,10mm



O.040
Jenkins GOUGE,
2,4,6,8,10mm



O.041
FARABEUF
Elevator Straight



O.041B
FARABEUF
Elevator Curved



O.042
Kerrison Punch
Upward Cutting



O.043
Lake Seeker



O.044
Houses Graft
Press Forceps



O.045
Lempert Mastoid Suction Tube
Size 1-2-3-4



O.046
Ferguson Suction
Tube Size 1-2-3-4



O.047
Verhoeven Suction Tube SWG 14,
15, 16 17, 18, 20, 21, 22, 23, 24




O.048
VERHOEVEN
Suction Tube
Holder



O.049
Tube Connector




O.050
Temporal Bone Holder




O.051



O.052




O.053




O.054

Micro Aural Crocodile
Forceps Serrated


- O.051 Fine Straight**
- O.052 Fine 1 & 2 Teeth**
- O.053 Right Curved**
- O.054 Left Curved**



O.055



O.056




O.057

Micro Aural Cup Forceps

- O.055 Straight Cup**
- O.056 Right Curved**
- O.057 Left Curved**



O.058



O.059



O.060

Micro Aural Scissors

- O.058 Right Curved**
- O.059 Left Curved**
- O.060 Straight Fine**



O.061
Mc-Gee Wire Closing
Forceps



O.062
Teflon Piston Holding
Forceps




O.063
Grommet Inserter



O.064
Piston Introducing Forceps



O.065 Up
O.066 Down



- O.067 WULLSTEIN Needle Straigh, Sharp
- O.068 WULLSTEIN Needle Slight Curved, Sharp
- O.069 WULLSTEIN Needle Strong Curved, Sharp
- O.070 WULLSTEIN Needle Very Slight Curved, Sharp
- O.071 MICRO PICK 90° Right Angled
- O.072 MICRO PICK 45° Angled
- O.073 PLESTER Foot Plate Hook
- O.074 SCHUHKNECHT Micro Needle Straigh, Sharp
- O.075 SCHUHKNECHT Micro Needle 90° Right Angled, Right
- O.076 SCHUHKNECHT Micro Needle 90° Right Angled, Left
- O.077 SCHUHKNECHT Micro Needle 45° Right Angled, Left
- O.078 SCHUHKNECHT Micro Needle 45° Right Angled, Right



- O.079 ROSEN Elevator Wide Length
- O.080 BEALES Raspatory
- O.081 Sickle Knife
- O.082 PLESTER First Incision Knife
- O.083 MYRINGOTOMY Knife



O.084
Rosen Curette



O.085
Rosen Curette
Double Curved



O.086
Rosen Curette
Curved



O.087
Circular Cutting
Knife 60°



O.088
Circular Cutting
Knife Curved



O.089
House Curette
Double Ended



O.090
LARKIN Hand
Trepine



O.091
LARKIN
Penetrator (Burr)



O.092
LARKIN
Perforator





O.093
House Measuring
Roads 3.5, 4.0,
4.5mm




O.094
Hook Ball
End Bent 45°



O.095
SHEA Piston
Depth Gauge
3.5, 4.0, 4.5mm




O.096
Fine Hook




O.097
Fine Knife
Downward
Cutting




O.098
Fine Knife
Upward
Cutting




O.099
Spear
Pointed
Fine Knife
Cutting Left



O.100
Spear
Pointed
Fine Knife
Cutting Right




O.101
Raspatory
Fine
Curved Left



O.102
Raspatory
Fine
Curved Right



O.103
Cutting Burr
SS



7.0



6.0



5.0



4.0



3.1



2.7



1.8

O.104
Tungstone
Carbide
Cutting Burr



7.0



6.0



5.0



4.0



3.1



2.7



1.8



O.105
Diamond
Burr



7.0



6.0



5.0



4.0



3.1



2.7



1.8




O.106
Oval Shape
Burr SS



7.0



6.0



5.0



BURR STAND BLOCK



O.108
Hand Piece
Straight




O.109
Hand Piece
Angled



O.110
Hand Piece
Spray



O.111
Hand Piece for



ENT SURGICAL MICROSCOPE (Make:GEM)

Serial No. O.113

Technical Specifications:-

- Head: Straight Binocular Tube
- Magnification: 0.4x, 0.6x, 1.0x, 1.6x, 2.5x
- Eye Pieces: SWF 12.5x paired super wide field
- Working Dist: F=200mm Objective Lens
- Fine Focusing: Manual Fine Focusing
- Light Source: 100W LED Light Source with bright white LED Light
- Arm: Counter balanced pantographic arm with 320° rotation
- Floor Stand: Mobile floor stand on four caster wheels for easy handling and absolute stability
- Power Supply: AC 220V – 240V

Optional Accessories:-

- Tilttable/Inclined Head
- Beam Splitter with 'C' Mount
- CCD Digital Camera
- LED Screen for Live Surgery
- Motorised with Foot Control
- Objective lens: 250mm, 300mm & 400mm



O.117 Sterilization and Storage Micro Instrument Case



O.118 Binocular Loupe



O.120 Micro Instrument Case with Silicone Mat



O.121 Rack for holding Aural Micro Instruments



O.122 Micro Instruments Box