



### Main Features:

#### Case:

Custom Designed Case  
Scratch Resistant and Solid Structure  
Lightweight

#### Endoscopy Camera:

*EndoView lite* Endoscopy Camera  
Sony CCD Image Sensor  
800tvl high resolution camera  
Anti-Glare Coated Lens

#### Light Source:

10W LED Light Source  
High Performance Bright LED  
Adjustable Brightness

#### Screen:

13.5" LED Panel  
HD Screen

### Specifications:

<b>Camera Resolution:</b>	<b>800tvl</b>
<b>Image Sensor:</b>	<b>1/3" Sony CCD Sensor</b>
<b>LED Power:</b>	<b>10W</b>
<b>LED Lux:</b>	<b>&gt;100000lux</b>
<b>Color Temperature:</b>	<b>6500k</b>
<b>Screen Size:</b>	<b>13.5"</b>
<b>Weight:</b>	<b>5kg.</b>