

# TEMPORAL BONE HOLDER

## MODELS: GM225 / GM325



# GEM SURGITECH

MANUFACTURER OF MEDICAL EQUIPMENTS & INSTRUMENTS



### Application:

An essential tool for temporal bone laboratories and teaching facilities.

### Main Features:

- Sharp Prongs Ensures Strong Grip
- Heavy Base Helps with Stability
- Stainless Steel Bowl in Black Finish
- 3 Adjustable Screw-In Clamps Securely Holds Temporal Bone for Dissection  
\*5 Adjustable Screw-In Clamps for model GM325\*
- Adjustable Bowl Angle and Lockable at Desired Positioning
- Outflow Funnel for Irrigation Liquid

### Dimensions:

Model:	GM225	GM325
Circumference of Inner Bowl:	535mm	580mm
Height of Bowl:	75mm	85mm
Total Height:	185mm	195mm